


Lightning Blot System




Move your research ahead faster with the Lightning Blot System. It transfers proteins from SDS-PAGE gels to membranes in 6-12 minutes, delivering quantitative results and improved reproducibility.

- **FAST:** 2-minute setup, complete electroblotting in <12 minutes
- **QUANTITATIVE TRANSFER** for proteins up to 150 kDa, good results to 250 kDa and higher
- **MORE REPRODUCIBLE** than standard tank transfers
- **GREEN:** Less packaging waste, no transfer buffer disposal
- **ON BUDGET:** Cost / blot similar to standard tank transfer

Quantitative Transfer of High and Low Molecular Weight Proteins



Chemiluminescence Detection of AKT in Cell Lysate after Lightning Blot Transfer



Lane 1: Protein Ladder

Lane 2-9: 2-fold dilution series of lysate prepared from HEK 293 cells using the AlphaScreen SureFire lysis buffer

Transfer performed using Lightning Blot semi-dry transfer system. Blocked for 1 hour with PerkinElmer blocking solution. Incubated with anti-AKT antibody at 1:1,000 dilution and anti-rabbit-HRP at 1:50,000 dilution. Chemiluminescence signal generated with **Western Lightning ECL Pro**.

	Lightning Blot System	Semi-dry Transfer	Tank Transfer
Buffer Preparation	0 minutes	15 minutes	15 minutes
Equilibration of Gel in Transfer Buffer	0 minutes	20 minutes	0 minutes
Transfer Sandwich Assembly	2 minutes	10 minutes	10 minutes
Transfer	6-12 minutes	30-90 minutes	1-3 hours
Cleanup	0 minutes	10 minutes	10 minutes
Total Transfer Time	8-14 minutes	1.5-2.5 hours	1.5-3.5 hours
Results	<ul style="list-style-type: none"> • No buffer required • Good transfers for proteins from 10 to 250 kDa • Quantitative transfer for proteins up to at least 150 kDa • Minimal "blow through" of low molecular weight proteins 	<ul style="list-style-type: none"> • Uses up to 250 mL of buffer • Often does not work well for high molecular weight proteins • Potential for short-circuit and failed transfer if membrane and filter paper do not match gel size 	<ul style="list-style-type: none"> • 1L of 20% methanol buffer for disposal after each run • Increasing temperature during a long run with high current may affect transfer efficiency and reproducibility

Ordering Information

Visit www.perkinelmer.com/lightningblot to order online.

Description	Dimensions	Catalog number
Lightning Blotter Mini Transfer System (for 1 miniblot)	10 cm x 10 cm	NEF2000001EA
Lightning Blotter Midi Transfer System (for 1 midiblot or 2 miniblots)	10 cm x 18 cm	NEF201001EA
Lightning Blot Mini Transfer Stack 2pk	8 cm x 8 cm	NEF211001EA
Lightning Blot Mini Transfer Stack 10pk	8 cm x 8 cm	NEF212001EA
Lightning Blot Mini Transfer Stack 50pk	8 cm x 8 cm	NEF213001EA
Lightning Blot Midi Transfer Stack 2pk	8.5 cm x 13.5 cm	NEF221001EA
Lightning Blot Midi Transfer Stack 10pk	8.5 cm x 13.5 cm	NEF222001EA
Lightning Blot Midi Transfer Stack 50pk	8.5 cm x 13.5 cm	NEF223001EA



PerkinElmer, Inc.
 940 Winter Street
 Waltham, MA 02451 USA
 P: (800) 762-4000 or
 (+1) 203-925-4602
www.perkinelmer.com



For a complete listing of our global offices, visit www.perkinelmer.com/ContactUs

Copyright ©2013 PerkinElmer, Inc. All rights reserved. PerkinElmer® is a registered trademark of PerkinElmer, Inc. All other trademarks are the property of their respective owners.

010968_03